

# Ambu<sup>®</sup> Mark III

## Main Features Include:

- Developed according to US military specifications
- Butyl outer cover for optimal resistance to penetration during ventilation in toxic atmospheres
- Can be used in all conventional ventilation emergencies
- NATIONAL STOCK NUMBER (NSN): #6515-01-338-6602
- PART NUMBER: #233-005-000



## Ambu<sup>®</sup> Mark III

Ambu is pleased to announce the availability of the world's only effective, truly portable manual resuscitation system designed and tested specifically for use in chemical or biological conditions.

### The Ambu rdic is supplied as a complete system including:

- A genuine Ambu Mark III reusable resuscitator which is equipped with impermeable elastomeric components on the compression bag and patient valve (butyl rubber)
- A highly flexible butyl rubber extension tube which offers unparalleled flexibility during casualty transport
- A contoured face mask with clear dome which permits better monitoring of the patient

- A specialized inlet filtration system, including C2 canister, for filtering contaminated air
- A durable, high-impact molded carry case



### Ambu S.r.l.

Via Paracelso, 20  
Centro Direzionale Colleoni  
Palazzo Andromeda 3  
20864 Agrate Brianza MB - IT  
T +39 039 657811  
F +39 039 6898177  
www.ambu.it

### Ambu A/S

Baltorpbakken 13  
DK - 2750 Ballerup  
Denmark  
T +45 72 25 20 00  
F +45 72 25 20 50  
www.ambu.com